


Derma Shine[®] Pure Mask

FACE LED MASK



Time for treatment :
9min, 4times continuous sue when fully charged.

- Yellow + IR : this mode helps restore facial skin by preventing the creation of fine wrinkles and thick wrinkles on the face.
- Red + IR : This mode stimulates the fibroblast activity of the facial skin and helps to regenerate collagen and elastin and to recover the wound quickly
- Green + IR : This mode induce cell growth and helps to repair and maintain skin, so it can relieve black spots on skin surface and lighten skin color



SPECIFICATION

Dimension (H X W) : 183 mm X 234 mm

Battery : 2,850 mAh

LED 408ea (LED With 4 different types of wavelength emission)

Green : 525 nm / 102 ea Red : 620 nm / 102 ea

Yellow : 590 nm / 102 ea IR (Infrared Ray) : 850 nm / 102 ea



(Sangeaewon-dong, Halla Sigma Valley) #405 #406,
545, Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13215, Korea
Tel. 82-70-7492-5955 Fax. 82-70-7469-5122